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

The Artificial Arch Foot System brings an exceptional level of comfort to the low to high activity user. Artificial arch spring enhances stability at heel strike through shock absorption and flexion as the artificial spring flexes and extends.

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Standard bench alignment techniques are used for initial alignment of the 1C70 Foot. The foot is designed with a 1 cm heel height. Align the center of the pyramid of the foot 10 mm anterior to the sagittal bisection of the proximal brim. Set the desired 5-7 degrees flexion angle in the sagittal plane. In the coronal plane set the desired 3-5 degrees adduction angle to accommodate for knee genu varum. The foot should also be set with 0-12 mm inset relative to the bisection of the proximal brim. Less inset is required for a shorter limb length. Rotate the foot with approximately 5 degrees of external rotation relative to the line of progression.

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\* If you are not sure which foot to select, please consult with one of our technical service representatives.

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USA YOBAND PROSTHETIC & ORTHOPAEDIC CO.,LTD